

BA Ultimate Hygiene Range

BA Ultimate Hygiene Range includes a wide range of products for all your needs: autoclavable treatment centres including handpieces and a range of tips covering scaling, periodontic and endodontic applications.



KEY FACTS:

BA Ultimate Autoclavable Treatment centre available for BA/EMS® type, Satelec® and Sirona® type handpieces. Each treatment centre contains a fully autoclavable handpiece, 2 scaling tip (BAC161, BAC164), 1 perio tip (BAC173), 1 endo head (BAC151), a tip changer and an endo wrench.

BA Ultimate Treatment Centre

Code	Description
BAC150E	Autoclavable treatment centre with BA/EMS® type handpiece
BAC150S	Autoclavable treatment centre with Satelec® type handpiece
BAC150A	Autoclavable treatment centre with Sirona® type handpiece
BAC150AL	Autoclavable treatment centre Sirona® type handpiece with light

BA Ultimate tip changer

Code	Description
BACS122	Tip changer for scaling & perio tips (EMS® Satelec®)
BACS127	Tip changer for scaling & perio tips (Sirona®)

BA Ultimate Tips (pack of 3 tips - please choose correct fitting*)

 <p>BAC161 Scaling Removes Supragingival Calculus</p>	 <p>BAC162 Scaling Removes pragingival Heavy Calculus NOTE: Not Available For Sirona Type Handpieces</p>	 <p>BAC163 Scaling Removes Calculus On Supragingival, Subgingival & Interdental NOTE: Not Available For Sirona Type Handpieces</p>	 <p>BAC164 Scaling Removes All Supragingival Calculus</p>	 <p>BAC165 Scaling Removes Supragingival & Neck Calculus</p>	 <p>BAC166 Scaling Removes Supragingival Heavy Calculus</p>
 <p>BAC171 Perio Removes Supragingival Calculus</p>	 <p>BAC173 Perio Subgingival perio use</p>	 <p>BAC174 Perio Removes Supragingival & Interdental Calculus</p>	 <p>BAC151 Endo 120° For Anterior Root Canal Cleaning</p>	 <p>BAC152 Endo 95° For Posterior Root Canal Cleaning</p>	 <p>BAC153 Endo For Lateral Condensation</p>

*Please suffix tip with 'E' for EMS and BA International type handpieces, 'S' for Satelec® and NSK® type handpieces, 'A' for Sirona® type handpieces e.g. BAC161E.