




NEW

VRN-Q6 ULTRASONIC + AIR POLISHER

ULTRASONIC PERIODONTAL THERAPY SYSTEM

700ml high capacity water bottle x 2 pcs, internal anti-corrosion liquid ways

New water regulating design, adjust the water at the handpiece.



New water regulating design, adjust the water at the handpiece.

3.5" color TFT LCD display

Waterproof IML touch operation panel.

Wireless foot control, exclusive patent.





Supragingival Air Polishing

Air polishing with sand blasting
Effectively remove the dental
plaque in Supragingival and
subgingival with sand blasting
which driven by high pressure
air to prevent the oral disease.

1



Subgingival Air Polishing

RESIDUAL POCKETS / DEEP TRANS
MUCOSAL PORTION
Preserves integrity of cementum,
dentin and implant surfaces
during treatment.

2



No pain Supra and
Subgingival Scaling

NO-PAIN ultrasonic piezo treatment
It combines scaling, periodontal and
endodontic in one device. 5 stainless
tips and 10 titanium alloy perio
treatment tips attached, they can
meet the most needs of supragingival
and subgingival scaling operations.

3



Implant Maintain

Dental implant maintaining
Sand blasting could remove the
plaque and stains on the implant
without any harm to the material.
Besides, 1 PEEK plastic special
implant maintaining tip could
offer the soft and high efficient
operation.

4